

SYPRO® Ruby Staining of ONCYTE® Porous Nitrocellulose Film Slides

1. Purpose

- 1.1. General staining of proteins deposited on porous nitrocellulose films for sample protein quantification.

2. Equipment & Materials

2.1. Equipment

- 2.1.1. Orbital Shaker
- 2.1.2. Microarray Fluorescence Scanner such as the GenePix 4000B.
Microarray Scanner (Molecular Devices)

2.2. Materials

- 2.2.1. Glass staining dishes with lids (Wheaton Cat# 900203).
- 2.2.2. Glass slide racks with handles (Wheaton Cat# 900204, 900205).
- 2.2.3. ProPlate™ modules (Grace Bio-Labs)

2.3. Reagents

- 2.3.1. SYPRO® Ruby (Invitrogen Cat# S12000)
- 2.3.2. Super G Blocking Reagent, optional (Grace Bio-Labs)
- 2.3.3. Fixative Solution: 7% acetic acid, 10% methanol in distilled water
- 2.3.4. Distilled water

3. Procedure

- 3.1. Optional: Submerge spotted protein array in Super G blocking reagent for 1 hour with no agitation.
 - 3.1.1. Use of Super G blocking reagent will provide lower background fluorescence levels.
- 3.2. Transfer slides to a staining jar containing fresh water.
 - 3.2.1. Wash 3 times for 15 minutes each, with agitation (105 rpm on orbital shaker).
- 3.3. Transfer slides to a staining jar containing the fixative solution.
 - 3.3.1. Incubate for 15 minutes, with agitation (105 rpm on orbital shaker).
- 3.4. Transfer slides to a staining jar containing fresh water.
 - 3.4.1. Wash 4 times for 5 minutes each, with agitation (105 rpm on orbital shaker).
- 3.5. Transfer slides to a staining jar containing the SYPRO® Ruby stain.
 - 3.5.1. Alternatively: To save on reagent, instead of a staining jar use a ProPlate™ module such as the 1-well module.
 - 3.5.2. Apply a 1-well ProPlate™ module to the wet slide (work quickly to prevent slide drying) and add 2 ml of SYPRO® Ruby stain.
- 3.6. Perform staining for 30 minutes at room temperature, with agitation (105 rpm on orbital shaker).
 - 3.6.1. If using a ProPlate™ module, perform incubation with agitation at 50 rpm.

SYPRO® Ruby Staining of ONCYTE® Porous Nitrocellulose Film Slides

- 3.7. Transfer slides to a staining jar containing fresh water.
 - 3.7.1. Wash 4 times for 1 minute each, with agitation (105 rpm on orbital shaker).
 - 3.7.2. If using a ProPlate™ module, remove SYPRO® Ruby from chamber with a pipette and replace with 2 ml fresh water. Perform the first water wash in the ProPlate™ module, then unclip and transfer slide to a staining jar for washes 2-4.
- 3.8. Air-dry the slides in the dark.
- 3.9. Scan slides using a fluorescent scanner (these settings serve as a starting point and may vary from scanner to scanner):
 - 3.9.1. Axon 4000B Scanner
 - 3.9.1.1. 532 nm: 33% laser power, 400 PMT setting
 - 3.9.1.2. 635 nm: 100% laser power, 800 PMT setting

See Appendix for typical results obtained with this method

4. References

1. SYPRO® Ruby Protein Gel Stain Product Information, Molecular Probes, Rev. 11/20/2007.
2. SYPRO® Ruby Protein Stains Instruction Manual, BioRad, 4006173 Rev B.

SYPRO® Ruby Staining of ONCYTE® Porous Nitrocellulose Film Slides

Appendix

Figure 1. SYPRO® Ruby Protein Quantification on ONCYTE® Film Slides.




Figure 2. Fluorescent Signal from SYPRO® Ruby Stained Protein on ONCYTE® Film Slides.




Figure 2. Fluorescent signal after staining with SYPRO® Ruby is linear over 3 orders of magnitude. Shown is the dynamic range of protein staining with SYPRO® Ruby on a SuperNOVA nitrocellulose film slide blocked with Super G blocking reagent. BSA standard was spotted in various RPPA lysate buffers with R^2 averaging 0.9955 for all of the buffers (Buffer A: SDS/Tris, Buffer B: SDS/Tris/Triton; Buffer C: Urea; Buffer D: Grace Bio-Labs RPPA lysis buffer).