


CoverWell™ Perfusion Chambers


Form removable, water-tight multiwell chambers by pressing to coverglass or microscope slides.
Thin (0.18 mm / 0.007") polycarbonate cover has excellent transparency.

622101-501




Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
PC1R-0.5	0.6	180	40
PC1R-1.0	0.9	300	40
PC1R-2.0	1.7	520	20
PC1R-2.5	2.5	750	20

622102-502




Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
PC1L-0.5	0.6	350	40
PC1L-1.0	0.9	550	40
PC1L-2.0	1.7	1000	20
PC1L-2.5	2.5	1500	20

622103-503



Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
PC3L-0.5	0.6	75	40
PC3L-1.0	0.9	120	40
PC3L-2.0	1.7	200	20
PC3L-2.5	2.5	300	20


622104-504



PC4L-

Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
PC4L-0.5	0.6	70	40
PC4L-1.0	0.9	110	40
PC4L-2.0	1.7	200	20
PC4L-2.5	2.5	300	20


622105-505



PC8R-

Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
PC8R-0.5	0.6	35	20
PC8R-1.0	0.9	60	20
PC8R-2.0	1.7	100	10
PC8R-2.5	2.5	150	10

622551



PC1S-

Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
622551	0.6	180	40