

CultureWell™ Silicone Gaskets



Isolate cells using removable hydrophobic barriers on coverglass.


Reusable CultureWell gaskets are available bulk packaged.

Gaskets are non-sterile and may be sterilized by autoclaving, UV or alcohol.

103210/103220

NON-STERILE

CW 2R-



CW 2R-1.0 ← 1.0 mm silicone


CW 2R-2.0 ← 1.5 mm silicone

Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
103210	1.0	250 - 400	10
103220	1.5	300 - 500	10

103230

NON-STERILE

CW 3S-1.0




Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
103230	1.0	300 - 500	10

103240

NON-STERILE

CW 4R-1.0




Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
103240	1.0	50 - 100	10

103280

NON-STERILE

CW 8R-1.0




Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
103280	1.0	15 - 30	10

103250

NON-STERILE

CW 50R-1.0



Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
103250	1.0	3 - 10	10