

Isolate cells using removable hydrophobic barriers on coverglass.



Each silicone gasket is pre-applied to a No. 1.5 desag coverglass. There are 4 gasketed coverglass per sterile pouch. There are 5 pouches per package for a total of 20 parts per order.

CultureWell silicone gaskets and coverglass (non-sterile) may be purchased separately - see below.

103310/103320

STERILE

CWCS 2R-



Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
----------	--------------------------	----------------------	-------------------



103310	1.0	250 - 400	20
--------	-----	-----------	----

103220	1.5	300 - 500	20
--------	-----	-----------	----

103330

STERILE

CWCS 3S-1.0





Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
----------	--------------------------	----------------------	-------------------

103330	1.0	300 - 500	20
--------	-----	-----------	----

103340

STERILE

CWCS 4R-1.0





Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
----------	--------------------------	----------------------	-------------------

103340	1.0	50 - 100	20
--------	-----	----------	----

103380

STERILE

CWCS 8R-1.0




Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
----------	--------------------------	----------------------	-------------------

103380	1.0	15 - 30	20
--------	-----	---------	----


103350

STERILE

CW 50R-1.0



CWCS 50R-1.0




Item No.	Silicone Well Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
103350	1.0	3 - 10	20

162450

NON-STERILE

2450-1.5



No. 1.5  
Thickness  
Coverglass Only

Item No.	Quantity Per Pack
162450	1 oz