

CoverWell™ Imaging Chambers




Form removable, water-tight chambers to support thick and free-floating specimens for imaging.

Thin (0.18 mm / 0.007") polycarbonate cover has excellent transparency.

For additional sealing or permanent mounting in aqueous media, select "Adhesive Imaging Chambers"

635011-051




PCI-

	Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
PCI-0.5	635051	0.6	180	40
PCI-1.0	635011	0.9	300	40
PCI-2.0	635021	1.7	520	20
PCI-2.5	635031	2.5	750	20

631011-051


Adhesive Imaging Chambers



PCI-A-

	Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
PCI-A-0.5	631051	0.8	180	40
PCI-A-1.0	631011	1.3	300	40
PCI-A-2.0	631021	2.3	520	20
PCI-A-2.5	631031	2.8	750	20

635041



PCI-PAF

Item No.	Chamber Depth (mm)	Volume Per Well (ul)	Quantity Per Pack
635041	0.26	180	25

