


## | CultureWell™ – Removable Chamber Slides, Sterile |

### Chamber Slides

Multi-well cell culture chamber slide lets the user grow cells directly onto a standard borosilicate glass slide with an easily removable polystyrene chamber. The chamber may be detached from the slide by simply disassembling the well attachment.



- Polypropylene bracket and lid, borosilicate glass slide and polystyrene chambers
- Gamma irradiation sterilization
- Easy disassembly for microscope viewing
- Non-pyrogenic
- Non-cytotoxic
- DNase and RNase free
- Available in 1, 2, 4, and 8 well configurations, see table below for details

| CultureWell™ – Removable Chamber Slides, Sterile |

Catalog Item #	Well	Chamber Color	Growth Area (cm <sup>2</sup> )	Well volume (mL)	Well dimension (mm)
103510	1	White	9.40	2.50 - 5.50	23.00 x 49.60
103511	1	Clear	9.40	2.50 - 5.50	23.00 x 49.60
103512	1	Black	9.40	2.50 - 5.50	23.00 x 49.60
103520	2	White	4.55	1.20 - 2.50	23.00 x 22.20
103521	2	Clear	4.55	1.20 - 2.50	23.00 x 22.20
103522	2	Black	4.55	1.20 - 2.50	23.00 x 22.20
103530	4	White	2.13	0.50 - 1.30	23.00 x 11.40
103531	4	Clear	2.13	0.50 - 1.30	23.00 x 11.40
103532	4	Black	2.13	0.50 - 1.30	23.00 x 11.40
103540	8	White	0.98	0.20 - 0.60	11.50 x 11.40
103541	8	Clear	0.98	0.20 - 0.60	11.50 x 11.40
103542	8	Black	0.98	0.20 - 0.60	11.50 x 11.40